

ASCP-MLT^{Q&As}

MEDICAL LABORATORY TECHNICIAN - MLT(ASCP)

Pass ASCP ASCP-MLT Exam with 100% Guarantee

Free Download Real Questions & Answers **PDF** and **VCE** file from:

<https://www.leads4pass.com/ascp-mlt.html>

100% Passing Guarantee
100% Money Back Assurance

Following Questions and Answers are all new published by ASCP
Official Exam Center

-  **Instant Download** After Purchase
-  **100% Money Back** Guarantee
-  **365 Days** Free Update
-  **800,000+** Satisfied Customers



QUESTION 1

The part of the microscope that control the amount of light entering the specimen much like the iris of your eye controls light is called the iris _____.

- A. fine adjustment
- B. objective
- C. condenser
- D. diaphragm

Correct Answer: D

QUESTION 2

Measuring the turnaround times of troponin tests from the emergency department would provide valuable information. This data represents circumstances in which patient safety may be compromised if results are delayed and treatment is not started as soon as possible.

A laboratory manager wants to evaluate the timeliness of patient services in order to prevent medical errors caused by delay in treatment. Measuring the turnaround time for which of these analytes would provide the most valuable information?

- A. Routine glucose tests from patients on medical floors
- B. Lipid panels from the outpatient clinic
- C. Troponin tests from the emergency department

Correct Answer: C

QUESTION 3

Match each of the descriptions with the appropriate magnification:

1.

Color, Rouleau, Overall Slide Quality, Cell Distribution

2.

Platelet estimates RBC-platelet-WBC morphology WBC differential RBC inclusions

3.

Select area to examine, WBC estimate

- A. 40X (Dry)

B. 10X

C. 100X (Oil)

Correct Answer: ABC

QUESTION 4

"Universal donor", (a misnomer) is usually applied to group O, Rh negative blood. Although it may be necessary to use group O, Rh negative blood in an extreme emergency, it is preferable to use type specific blood for emergencies. In an extreme emergency, if the ABO and Rh type are unknown which of the following should be given to the patient?

A. Group O, Rh positive blood

B. Group AB, Rh negative blood

C. Group O, Rh negative blood

D. Any blood type is OK

Correct Answer: C

QUESTION 5

Anti-A, anti-P, anti-Leb, and anti-M all react best at 40 C as they are predominantly IgM antibodies. Other antibody group choices above include IgG antibodies such as anti-K, anti-s, anti-S, and anti-Fya, anti-Lub, etc. which react best at 37o

C.

Which of the following groups of antibodies generally reacts most strongly at 40 C:

A. Anti-A, Anti-P1, Anti-Leb, Anti-M

B. Anti-B, Anti-K, Anti-Lua, Anti-Fya

C. Anti-H, Anti-S, Anti-Jkb, Anti-Leb

D. Anti-A, Anti-K, Anti-Lub, Anti-s

Correct Answer: A

[Latest ASCP-MLT Dumps](#)

[ASCP-MLT PDF Dumps](#)

[ASCP-MLT Practice Test](#)